


GOAT ANTI-FOXP3 ANTIBODY

SKU: EB07034



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	OTTMUSP00000017961 OTTMUSP00000017952 RP23-54C14.1 scurfy sf JM2, scurfin MGI:1891436 forkhead box P3 [Mus musculus] Foxp3
Accession ID	NP_473380.1
Blocking Peptide	EBP07034
Immunogen	Peptide with sequence CTFPRSGTPRKDSN, from the internal region of the protein sequence according to NP_473380.1.
Peptide Sequence	CTFPRSGTPRKDSN
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated

Predicted Species	Mouse, Rat
Reactive Species	Human
Mouse Gene ID	20371
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Spleen. Recommended concentration: 5µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:16000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

