

# GOAT ANTI-FOXQ1 / FKHRL1 ANTIBODY

**SKU:** EB06032

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

## SPECIFICATIONS

**Formulation** Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 µg

**Storage Instructions** Aliquot and store at -20°C. Minimize freezing and thawing.

**Synonym /**

**Alias** winged helix/forkhead transcription factor|forkhead box Q1|HFX1|FOXQ1

**Names**

**Accession ID** NP\_150285.3

**Blocking Peptide** EBP06032

**Immunogen** Peptide with sequence KLEVFVPRAAHGD-C, from the N Terminus of the protein sequence according to NP\_150285.3.

**Peptide Sequence** KLEVFVPRAAHGD-C

**Purification Method** Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

**Shipping Instructions** Refrigerated

**Predicted Species** Human

**Reactive Species** Human

**Human Gene ID** 94234

**Mouse Gene ID** 15220

**Product Grade** [https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\\_medium.png](https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png)

**ELISA Detection Limit** Antibody detection limit dilution 1:32000.

**Western Blot** Approx 50kDa band (faint) observed in H460 lysate (predicted MW of 46kDa according to NP\_150285). Recommended for use at 1-3µg/ml.

**Application Type** Pep-ELISA, WB

## DOCUMENTS

- [Data Sheet](#)

## GALLERY IMAGES

