


GOAT ANTI-FPRL1 (AA328-339) ANTIBODY

SKU: EB06979



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	formyl peptide receptor 2 FPR2 lipoxin A4 receptor (formyl peptide receptor related) formyl peptide receptor-like 1 FMLP-R-II LXA4R FPRH2 FPRH1 FPR2A FMLPX HM63 ALXR FPRL1
Usage Summary	Flow Cytometry: Flow cytometric analysis of K562 cells. Recommended concentration: 10ug/ml.
Accession ID	NP_001005738.1; NP_001453.1
Blocking Peptide	EBP06979
Immunogen	Peptide with sequence C-DSAPTNDTAANS, from the C Terminus of the protein sequence according to NP_001005738.1; NP_001453.1.
Product Comments	Both reported variants (NP_001005738.1 and NP_001453.1) are identical proteins.

Peptide Sequence	C-DSAPTNDTAANS
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	2358
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:1000.
Application Type	Pep-ELISA, FC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

