

# GOAT ANTI-FRA1 / FOSL1 ANTIBODY

**SKU:** EB07279



250kDa  
150kDa  
100kDa  
75kDa  
50kDa  
37kDa  
25kDa  
20kDa  
15kDa  
10kDa

## SPECIFICATIONS

**Formulation** Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 µg

**Storage Instructions** Aliquot and store at -20°C. Minimize freezing and thawing.

**Synonym /**

**Alias** FOS-like antigen-1|fra-1|HGNC:13718|FOS-like antigen 1|FRA1|FOSL1

**Names**

**Accession ID** NP\_005429.1

**Blocking Peptide** EBP07279

**Immunogen** Peptide with sequence C-QREIEELQKQKER, from the internal region of the protein sequence according to NP\_005429.1.

**Product Comments** No cross-reactivity is expected with V-FOS.

**Peptide Sequence** C-QREIEELQKQKER

**Purification Method** Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

**Shipping Instructions** Refrigerated

**Predicted Species** Human, Mouse, Rat, Dog, Pig, Cow

**Reactive Species** Human

**Human Gene ID** 8061

**Mouse Gene ID** 14283

**Rat Gene ID** 25445

**Product Grade** [https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\\_medium.png](https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png)

**ELISA**

**Detection Limit** Antibody detection limit dilution 1:1000.

**Western Blot** Approx 26-28kDa band observed in Human Bone Marrow, Breast and Placenta lysates (calculated MW of 29.4kDa according to NP\_005429.1). Recommended concentration: 0.1-0.3µg/ml.

**Application Type** Pep-ELISA, WB

## DOCUMENTS

- [Data Sheet](#)

## GALLERY IMAGES

