

GOAT ANTI-FRAT2 / GSK-3 ANTIBODY

SKU: EB05377



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	FRAT2 GSK-3 frequently rearranged in advanced T-cell lymphomas 2 GSK-3 binding protein FRAT2 RP11-452K12.2 MGC10562
Accession ID	NP_036215.1
Blocking Peptide	EBP05377
Immunogen	Peptide with sequence C-GPDRIALQPSGSLL, from the C Terminus of the protein sequence according to NP_036215.1.
Peptide Sequence	C-GPDRIALQPSGSLL
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	23401
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Approx 30-35kDa band seen in Human Heart lysate. Recommended for use at 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

