Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-FRIZZLED 3 ANTIBODY

SKU: EB08430

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias frizzled 3|hFz3|Fz-3|frizzled homolog 3 (Drosophila)|FZD3

Names

Accession NP 059108.1 ID

Blocking

EBP08430 **Peptide**

Peptide with sequence C-DDQRSKAGSIHSK, from the internal region of the protein sequence according to Immunogen

NP 059108.1.

Peptide Sequence

C-DDQRSKAGSIHSK

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human, Mouse, Rat

Human

7976 Gene ID

Mouse

14365

Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

Grade **ELISA**

Detection

Antibody detection limit dilution 1:2000.

Limit

Blot

Western

Preliminary experiments in Human Brain (Frontal Cortex), Mouse Brain and Rat Brain lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from

people in the field - have any results been reported with other antibodies/lysates?

Application

Pep-ELISA **Type**





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600