

GOAT ANTI-FRIZZLED 4 ANTIBODY

SKU: EB08356



250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names frizzled 4|exudative vitreoretinopathy 1 (autosomal dominant Criswick-Schepens syndrome)|WNT receptor frizzled-4|MGC34390|GPCR|FzE4|Fz-4|FZD4S|FEVR|EVR1|FZD4|frizzled homolog 4 (Drosophila)|Frizzled 4

Accession ID NP_036325.2

Blocking Peptide EBP08356

Immunogen Peptide with sequence DEEERRCDPIRIS, from the internal region of the protein sequence according to NP_036325.2.

Peptide Sequence DEEERRCDPIRIS

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human, Rat

Reactive Species Human

Human Gene ID 8322

Rat Gene ID 64558

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA Detection Limit Antibody detection limit dilution 1:16000.

Western Blot Approx 60kDa band observed in Human Colon and Small Intestine (Ileum) lysates (calculated MW of 59.9kDa according to NP_036325.2). Recommended concentration: 0.3-1µg/ml.

Application Type Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

