

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-FRIZZLED 4 ANTIBODY

SKU: EB08356



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa 75kDa 50kDa 37kDa

25kDa 20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

frizzled 4|exudative vitreoretinopathy 1 (autosomal dominant Criswick-Schepens syndrome)|WNT receptor frizzled-4|MGC34390|GPCR|FzE4|Fz-4|FZD4S|FEVR|EVR1|FZD4|frizzled homolog 4 (Drosophila)|Frizzled 4

Accession

Names

NP_036325.2

Blocking

EBP08356 **Peptide**

Peptide with sequence DEEERRCDPIRIS, from the internal region of the protein sequence according to **Immunogen**

NP 036325.2.

Peptide

Sequence

DEEERRCDPIRIS

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Rat

Species Reactive

Species

Human

Gene ID

Human 8322

Rat Gene ID 64558

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. a mazon aws.com/assets/products/image/elite_medium.png$

Grade **ELISA**

Detection Antibody detection limit dilution 1:16000.

Limit

Western

Approx 60kDa band observed in Human Colon and Small Intestine (Ileum) lysates (calculated MW of 59.9kDa

according to NP_036325.2). Recommended concentration: 0.3-1µg/ml.

Application

Type

Blot

Pep-ELISA, WB

GALLERY IMAGES





Telephone: (650) 697-3600

250kDa 150kDa 100kDa 75kDa 50kDa

e erest

37kDa

25kDa

20kDa

15kDa