



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-FYB ANTIBODY

SKU: EB05205

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

FYN-binding protein (FYB-120/130)|SLAP-130|ADAP|FYN binding protein (FYB-120/130)|FYB Alias

Names

Accession NP 001456

Immunogen Peptide with sequence YDDIADGCIYDND, from the C Terminus of the protein sequence according to NP_001456.

Peptide YDDIADGCIYDND

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Instructions Refrigerated

Predicted

Human, Mouse

Species

Human

2533 **Gene ID**

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade

Western

Western blot: No signal obtained yet but low background observed in Jurkat and U937 lysates at upto 2µg/ml. Blot