



GOAT ANTI-FYB ANTIBODY

SKU: EB05205



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	FYN-binding protein (FYB-120/130) SLAP-130 ADAP FYN binding protein (FYB-120/130) FYB
Accession ID	NP_001456
Immunogen	Peptide with sequence YDDIADGCIYDND, from the C Terminus of the protein sequence according to NP_001456.
Peptide Sequence	YDDIADGCIYDND
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Human Gene ID	2533
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
Western Blot	Western blot: No signal obtained yet but low background observed in Jurkat and U937 lysates at upto 2µg/ml.