






Telephone: (650) 697-3600

## **GOAT ANTI-FZD9 ANTIBODY**

**SKU:** EB08508



## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

**Alias** 

frizzled 9|FZD3|frizzled homolog 9 (Drosophila)|FZD9

**Names** 

**Accession** 

**Blocking** 

ID

**Peptide** 

EBP08508

NP\_003499.1

Peptide with sequence C-TKTDPSLENPTH, from the C Terminus of the protein sequence according to

Immunogen

NP\_003499.1.

**Peptide** Sequence

C-TKTDPSLENPTH

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions





Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

**Predicted** 

Human, Mouse, Rat

**Species** Reactive

**Species** 

Human

Human Gene ID

8326

**Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

 $\textbf{IHC Results} \ \ \text{Paraffin embedded Human Testis and Thyroid. Recommended concentration: } 5 \mu g/ml.$ 

**ELISA** 

Detection

Antibody detection limit dilution 1:32000.

Limit

Application

Type

Pep-ELISA, IHC

## **GALLERY IMAGES**

