


GOAT ANTI-GABBR1 ANTIBODY

SKU: EB09093



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage	
Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	seven transmembrane helix receptor gamma-aminobutyric acid (GABA) B receptor 1 OTTHUMP00000109099 OTTHUMP00000029116 OTTHUMP00000029115 OTTHUMP00000029114 GABAB, subunit 1c GABA-B receptor hGB1a dJ271M21.1.2 dJ271M21.1.1 GPRC3A GABBR1-3 GABABR1 GABAB(1e) FLJ92613 gamma-aminobutyric acid (GABA) B receptor, 1 GABBR1
Accession ID	NP_001461.1; NP_068703.1; NP_068704.2; NP_068705.2
Blocking Peptide	EBP09093
Immunogen	Peptide with sequence DRLSCDGSRVHLLYK, from the C Terminus of the protein sequence according to NP_001461.1; NP_068703.1; NP_068704.2; NP_068705.2.
Product Comments	This antibody is expected to recognize reported isoforms a, b and c (NP_001461.1, NP_068703.1, NP_068704.2 resp.).
Peptide Sequence	DRLSCDGSRVHLLYK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Reactive Species	Human
Human Gene ID	2550
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	In paraffin embedded Human Cortex shows axonal staining. Recommended concentration: 5-8µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:4000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

