

GOAT ANTI-GABRA6 ANTIBODY

SKU: EB09103



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	gamma-aminobutyric acid A receptor, alpha 6 GABA subunit A receptor alpha 6 MGC116904 MGC116903 gamma-aminobutyric acid (GABA) A receptor, alpha 6 GABRA6
Accession ID	NP_000802.2
Blocking Peptide	EBP09103
Immunogen	Peptide with sequence C-HPDSKYHLKKR, from the internal region of the protein sequence according to NP_000802.2.
Peptide Sequence	C-HPDSKYHLKKR
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Reactive Species	Human
Human Gene ID	2559
Mouse Gene ID	14399
Rat Gene ID	29708
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	Approx. 55kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 51.0kDa according to NP_000802.2). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

GALLERY IMAGES

