


GOAT ANTI-GALNS ANTIBODY

SKU: EB08057



250kDa
150kDa
100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	chondroitinase N-acetylgalactosamine-6-sulfatase MPS4A GAS GALNAC6S galactosamine (N-acetyl)-6-sulfate sulfatase (Morquio syndrome, mucopolysaccharidosis type IVA) GALNS
Accession ID	NP_000503.1
Blocking Peptide	EBP08057
Immunogen	Peptide with sequence C-TTHNLEDHTKLP, from the internal region of the protein sequence according to NP_000503.1.
Peptide Sequence	C-TTHNLEDHTKLP
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	2588
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Approx 60kDa band observed in Human Bone Marrow and Placenta lysates (calculated MW of 58.0kDa according to NP_000503.1). Recommended concentration: 0.01-0.03µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

