



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-GALNS ANTIBODY

SKU: EB08057



Telephone: (650) 697-3600



250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias Names

chondroitinase|N-acetylgalactosamine-6-sulfatase|MPS4A|GAS|GALNAC6S|galactosamine (N-acetyl)-6-sulfate

sulfatase (Morquio syndrome, mucopolysaccharidosis type IVA)|GALNS

Accession

ID

NP_000503.1

Blocking

EBP08057 **Peptide**

Peptide with sequence C-TTHNLEDHTKLP, from the internal region of the protein sequence according to **Immunogen**

NP 000503.1.

Peptide

Sequence

C-TTHNLEDHTKLP

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human

Human

2588

Reactive

Species

Human

Gene ID

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection

Antibody detection limit dilution 1:64000.

Limit

Western **Blot**

Approx 60kDa band observed in Human Bone Marrow and Placenta lysates (calculated MW of 58.0kDa

according to NP_000503.1). Recommended concentration: 0.01-0.03µg/ml.

Application

Type

Pep-ELISA, WB

GALLERY IMAGES



Telephone: (650) 697-3600



250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa