

GOAT ANTI-GALP ANTIBODY

SKU: EB07756



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	GALP galanin-like peptide precursor galanin-like peptide, alarin
Accession ID	NP_149097.1
Blocking Peptide	EBP07756
Immunogen	Peptide with sequence C-HLPQMGDQDGKRE, from the internal region of the protein sequence according to NP_149097.1.
Peptide Sequence	C-HLPQMGDQDGKRE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	85569
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Preliminary experiments showed an approx. 20+28kDa band in HeLa cell lysates at 1ug/ml (calculated MW of 12.5kDa according to Human NP_149097.1).
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)