





Telephone: (650) 697-3600

GOAT ANTI-GALP ANTIBODY

SKU: EB07756

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias GALP|galanin-like peptide precursor|galanin-like peptide, alarin

Names

Accession NP 149097.1 ID

Blocking

EBP07756 **Peptide**

Peptide with sequence C-HLPQMGDQDGKRE, from the internal region of the protein sequence according to Immunogen

NP 149097.1.

Peptide C-HLPQMGDQDGKRE Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Human

85569 Gene ID

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

ELISA

Detection Antibody detection limit dilution 1:64000.

Limit

Western Preliminary experiments showed an approx. 20+28kDa band in HeLa cell lysates at 1ug/ml (calculated MW of

Blot 12.5kDa according to Human NP_149097.1).

Application

Pep-ELISA **Type**