



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-GALR2 ANTIBODY

SKU: EB07564



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa 20kDa

15kDa

10kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias GALR2|galanin receptor 2|GALNR2|MGC125983|MGC125984

Names

Accession

NP_003848.1 ID

Blocking

EBP07564 **Peptide**

Peptide with sequence C-RYPLHSRELRTPRN, from the internal region of the protein sequence according to **Immunogen**

NP 003848.1.

Peptide

C-RYPLHSRELRTPRN Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Rat, Dog, Cow **Species**

Reactive

Human **Species**

Human

8811 Gene ID

Mouse

14428 **Gene ID**

Rat Gene ID 29234

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade ELISA

Detection Antibody detection limit dilution 1:4000.

Limit

Western Approx 40kDa band observed in Human Heart lysates (calculated MW of 41.7kDa according to NP 003848.1).

Blot Recommended concentration: 0.3-1µg/ml.

Application

Pep-ELISA, WB

Type



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GALLERY IMAGES

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa