

GOAT ANTI-GARNL4 ANTIBODY

SKU: EB05024



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias	KIAA1039 protein GTPase activating RANGAP domain-like 4
Names	
Immunogen	Peptide with sequence C-FDKLSHASSGAGH, from the C Terminus of the protein sequence according to .
Peptide Sequence	C-FDKLSHASSGAGH
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	23108
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Western Blot: No signal obtained yet but low background observed in Jurkat and Brain lysates at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?