

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-GARNL4 ANTIBODY

SKU: EB05024

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias KIAA1039 protein|GTPase activating RANGAP domain-like 4

Names

Immunogen Peptide with sequence C-FDKLSHASSGAGH, from the C Terminus of the protein sequence according to .

Peptide

C-FDKLSHASSGAGH Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Instructions Refrigerated

Predicted

Human **Species**

Human

23108 **Gene ID**

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

Grade **ELISA**

Detection

Limit

Antibody detection limit dilution 1:2000.

Western

Blot

Western Blot: No signal obtained yet but low background observed in Jurkat and Brain lysates at up to 1µg/ml. We

would appreciate any feedback from people in the field - have any results been reported with other

antibodies/lysates?