



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-GATA3 ANTIBODY


SKU: EB05570

Telephone: (650) 697-3600















SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

enhancer-binding protein GATA3|GATA-binding protein 3|MGC5445|MGC5199|MGC2346|HDR|GATA binding **Alias**

protein 3|GATA3 Names

Accession

NP_001002295.1; NP_002042.1 ID

Blocking

EBP05570 **Peptide**

Peptide with sequence EVTADQPRWVSHH-C, from the N Terminus of the protein sequence according to **Immunogen**

NP_001002295.1; NP_002042.1.

Product Comments

This antibody is expected to recognize both reported isoforms (NP 001002295.1 and NP 002042.1).

Peptide

EVTADQPRWVSHH-C Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Rat, Dog, Pig, Cow **Species**

Reactive

Human, Mouse **Species**

Human

2625 Gene ID

Mouse

14462 **Gene ID**

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade FLISA

Antibody detection limit dilution 1:64000. Detection

Limit

Blot

Western Approx 46kDa band observed in Human Tonsil and in Mouse Kidney lysates. (calculated MW of 47.9kDa

according to NP_002042.1 and Mouse NP_032117.1). Recommended concentration: 0.3-1µg/ml.

Application

Pep-ELISA, WB **Type**





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GALLERY IMAGES

