

GOAT ANTI-GDF5 ANTIBODY

SKU: EB09364



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	cartilage-derived morphogenetic protein-1 SYNS2 LAP4 CDMP1 growth differentiation factor 5 GDF5
Accession ID	NP_000548.1
Blocking Peptide	EBP09364
Immunogen	Peptide with sequence C-NEIKARSGQDDKT, from the internal region of the protein sequence according to NP_000548.1.
Peptide Sequence	C-NEIKARSGQDDKT
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Horse
Human Gene ID	8200
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)