



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-GDF5 ANTIBODY**

**SKU:** EB09364

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

cartilage-derived morphogenetic protein-1|SYNS2|LAP4|CDMP1|growth differentiation factor 5|GDF5 Alias

**Names** 

Accession NP 000548.1 ID

**Blocking** 

EBP09364 **Peptide** 

**Immunogen** NP\_000548.1. Peptide with sequence C-NEIKARSGQDDKT, from the internal region of the protein sequence according to

**Peptide** C-NEIKARSGQDDKT Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Human, Mouse, Horse **Species** 

Human

8200 **Gene ID** 

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Grade **ELISA** 

**Detection** Antibody detection limit dilution 1:128000.

Limit

**Application** Pep-ELISA **Type**