


GOAT ANTI-GGT1 ANTIBODY

SKU: EB11992



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	GGT1 gamma-glutamyltransferase 1 CD224 D22S672 D22S732 GGT GGT 1 GTG gamma-glutamyl transpeptidase gamma-glutamyltranspeptidase 1 glutamyl transpeptidase glutathione hydrolase 1
Accession ID	NP_005256.2
Blocking Peptide	EBP11992
Immunogen	Peptide with sequence C-PSASKEPDNHVYTR, from the internal region (near N terminus) of the protein sequence according to NP_005256.2.
Product Comments	Reported variants represent identical protein: NP_038347.2, NP_001027536.1, NP_001027537.1, NP_005256.2.
Peptide Sequence	C-PSASKEPDNHVYTR
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	2678
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 60kDa band observed in Human Kidney lysates (calculated MW of 61.4kDa according to NP_005256.2). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES



