


GOAT ANTI-MICS1 / GHITM (AA117-127) ANTIBODY

SKU: EB09816



250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	dermal papilla-derived protein 2 DERP2 DKFZp566C0746 FLJ26584 growth hormone inducible transmembrane protein HSPC282 MICS1 mitochondrial morphology and cristae structure 1 My021 PTD010 TMBIM5 transmembrane BAX inhibitor motif containing GHITM
Accession ID	NP_055209.2
Blocking Peptide	EBP09816
Immunogen	Peptide with sequence C-PQYVKDRIHST, from the internal region of the protein sequence according to NP_055209.2.
Peptide Sequence	C-PQYVKDRIHST
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Rat, Dog, Cow
Reactive Species	Human
Human Gene ID	27069
Rat Gene ID	290596
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Western Blot	Approx 38kDa band observed in Human Brain (Cerebral Cortex and Cerebellum) lysates (calculated MW of 37.2kDa according to NP_055209.2). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

