

GOAT ANTI-GIPC3 ANTIBODY

SKU: EB06195



250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names FLJ40925|DKFZp686j1198|C19orf64|GIPC PDZ domain containing family, member 3|PDZ domain protein GIPC3|GIPC3

Accession ID NP_573568.1

Blocking Peptide EBP06195

Immunogen Peptide with sequence EGAAAREARGTET-C, from the N Terminus of the protein sequence according to NP_573568.1.

Peptide Sequence EGAAAREARGTET-C

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human

Reactive Species Human

Human Gene ID 126326

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA Detection Limit Antibody detection limit dilution 1:32000.

Western Blot Approx 40kDa band observed in human lymph node lysates (predicted MW of 36kDa according to NP_573568). Recommended for use at 0.2-1.5µg/ml.

Application Type Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

