

# GOAT ANTI-GIPR ANTIBODY

**SKU:** EB09744



---

## SPECIFICATIONS

<b>Unit Size</b>	100 µg
<b>Synonym / Alias Names</b>	MGC126722 HGNC:4271 gastric inhibitory polypeptide receptor GIPR
<b>Accession ID</b>	NP_000155.1; NP_001295347.1
<b>Blocking Peptide</b>	EBP09744
<b>Immunogen</b>	Peptide with sequence C-EEQRQLPERAFR, from the internal region (near C-Terminus) of the protein sequence according to NP_000155.1; NP_001295347.1.
<b>Peptide Sequence</b>	C-EEQRQLPERAFR
<b>Shipping Instructions</b>	Refrigerated
<b>Predicted Species</b>	Human
<b>Human Gene ID</b>	2696
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png</a>
<b>ELISA Detection Limit</b>	Antibody detection limit dilution 1:128000.
<b>Application Type</b>	Pep-ELISA

## DOCUMENTS

- [Data Sheet](#)