



Telephone: (650) 697-3600



GOAT ANTI-GIPR ANTIBODY

SKU: EB09744

SPECIFICATIONS

Unit Size 100 μg

Synonym /

Alias MGC126722|HGNC:4271|gastric inhibitory polypeptide receptor|GIPR

Names

ID

Accession

NP_000155.1; NP_001295347.1

Blocking

EBP09744 **Peptide**

Peptide with sequence C-EEQRQLPERAFR, from the internal region (near C-Terminus) of the protein sequence

according to NP_000155.1; NP_001295347.1.

Peptide

Sequence

C-EEQRQLPERAFR

Shipping

Instructions Refrigerated

Predicted

Species

Human

Human

Gene ID

2696

Product Grade

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. a mazon aws.com/assets/products/image/aspiring_medium.png$

ELISA

Detection Antibody detection limit dilution 1:128000.

Application Pep-ELISA **Type**