
Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-GLI1 ANTIBODY

SKU: EB07027



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias Names glioma-associated oncogene homolog 1 (zinc finger protein)|glioma-associated oncogene family zinc finger 1|GLI family zinc finger 1|Tax Helper Protein|glioma-associated oncogene homolog 1|glioma-associated oncogene homolog (zinc finger protein)|GLI|HGNC:4317|glioma-associated oncogene homolog 1 (zinc finger

protein)|GLI1

Usage

Immunofluorescence: Strong expression of the protein seen in the cytoplasm and nuclei of

U2OS cells. Recommended concentration: 10µg/ml. **Summary**

Accession

NP_005260.1; NP_001153517.1; NP_001161081.1

Blocking Peptide

EBP07027









Immunogen

Peptide with sequence C-EDAPAQEPSYQSPK, from the internal region of the protein sequence according to NP_005260.1; NP_001153517.1; NP_001161081.1.

Peptide

C-EDAPAQEPSYQSPK

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method

using the immunizing peptide.

Shipping Instructions

Refrigerated

Predicted

Species

Reactive

Species

Human

Human

Human

2735 Gene ID

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection

Antibody detection limit dilution 1:8000.

Limit

Application

Type

Pep-ELISA, IF

GALLERY IMAGES

