


Telephone: (650) 697-3600

GOAT ANTI-GLUCOGON RECEPTOR ANTIBODY





SPECIFICATIONS

 $\textbf{Unit Size} \qquad 100 \; \mu g$

Synonym /

Alias MGC138246|GGR|glucagon receptor|GCGR

Names

Flow cytometric analysis of HeLa cells. Recommended concentration: 10ug/ml.

Accession NP_

Blocking

NP_000151.1

Peptide

EBP07834

Immunogen

Peptide with sequence C-HGPPSKELQFGR, from the internal region (near the C Terminus) of the protein sequence according to NP_000151.1.

Peptide

Sequence

C-HGPPSKELQFGR





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Shipping Refrigerated Instructions

Predicted

Human **Species**

Reactive

Human

Species Human

2642

Gene ID Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws.com/assets/products/image/elite_medium.png$

Grade **ELISA**

Detection Antibody detection limit dilution 1:32000.

Limit

Type

 $\begin{array}{c} \textbf{Application} \\ \textbf{Pep-ELISA, FC} \end{array}$

