


GOAT ANTI-GLUCOGON RECEPTOR ANTIBODY

SKU: EB07834



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias / Names	MGC138246 GGR glucagon receptor GCGR
Usage Summary	<p>Immunofluorescence: A customer reported staining in the endoplasmic reticulum, cytoplasm and nucleus of U2OS cells at a concentration of 10ug/ml.</p> <p>Flow Cytometry: Flow cytometric analysis of HeLa cells. Recommended concentration: 10ug/ml.</p>
Accession ID	NP_000151.1
Blocking Peptide	EBP07834
Immunogen	Peptide with sequence C-HGPPSKELQFGR, from the internal region (near the C Terminus) of the protein sequence according to NP_000151.1.
Peptide Sequence	C-HGPPSKELQFGR

Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	2642
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA, FC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

