

GOAT ANTI-GNRPX / PLEKHJ1 ANTIBODY

SKU: EB06141



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	guanine nucleotide releasing protein x 9530063M10Rik likely ortholog of mouse guanine nucleotide releasing protein x FLJ10297 GNRPX pleckstrin homology domain containing, family J member 1 PLEKHJ1
Accession ID	NP_060519.1
Blocking Peptide	EBP06141
Immunogen	Peptide with sequence C-SEEARFQLSGLQA, from the C Terminus of the protein sequence according to NP_060519.1.
Peptide Sequence	C-SEEARFQLSGLQA
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	55111
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)