

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-GNRPX / PLEKHJ1 ANTIBODY

SKU: EB06141

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

guanine nucleotide releasing protein x|9530063M10Rik|likely ortholog of mouse guanine nucleotide releasing **Alias** protein x|FLJ10297|GNRPX|pleckstrin homology domain containing, family J member 1|PLEKHJ1

Names

Accession NP 060519.1 ID

Blocking

EBP06141 **Peptide**

Peptide with sequence C-SEEARFQLSGLQA, from the C Terminus of the protein sequence according to Immunogen

NP 060519.1.

Peptide

C-SEEARFQLSGLQA Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted Human

Species Human

Gene ID

55111

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

ELISA

Detection Antibody detection limit dilution 1:16000.

Limit

Application Pep-ELISA **Type**