


GOAT ANTI-GOLGA4 / GOLGIN-245 ANTIBODY

SKU: EB06924



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names tGolgin-1|trans-Golgi p230|golgin-245|golgin-240|72.1 protein|256 kDa golgin|p230|MU-RMS-40.18|GOLG|GCP2|HGNC:4427|golgi autoantigen, golgin subfamily a, 4|GOLGA4

Accession ID NP_002069.2

Blocking Peptide EBP06924

Immunogen Peptide with sequence C-DDQTQKILERED, from the C Terminus of the protein sequence according to NP_002069.2.

Peptide Sequence C-DDQTQKILERED

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species	Human, Mouse, Rat
Reactive Species	Human
Human Gene ID	2803
Mouse Gene ID	54214
Rat Gene ID	501069
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	In paraffin embedded Human Liver shows textured staining of the cytoplasm consistent with Golgi structures in hepatocytes. Recommended concentration, 3-6µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:4000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

