

# GOAT ANTI-GPR39 ANTIBODY

SKU: EB07567



## SPECIFICATIONS

<b>Unit Size</b>	100 µg
<b>Synonym / Alias Names</b>	G protein-coupled receptor 39 GPR39
<b>Accession ID</b>	NP_001499.1
<b>Blocking Peptide</b>	EBP07567
<b>Immunogen</b>	Peptide with sequence C-QSEAEPQSKSQ, from the internal region of the protein sequence according to NP_001499.1.
<b>Peptide Sequence</b>	C-QSEAEPQSKSQ
<b>Shipping Instructions</b>	Refrigerated
<b>Predicted Species</b>	Human
<b>Human Gene ID</b>	2863
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png</a>
<b>ELISA Detection Limit</b>	Antibody detection limit dilution 1:1000.
<b>Western Blot</b>	Preliminary experiments in Human Brain (cerebellum) lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
<b>Application Type</b>	Pep-ELISA

## DOCUMENTS

- [Data Sheet](#)

