

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-GPR39 ANTIBODY

SKU: EB07567

SPECIFICATIONS

Unit Size 100 μg

Synonym /

Alias G protein-coupled receptor 39|GPR39

Names

Accession NP_001499.1 ID

Blocking

EBP07567 **Peptide**

Immunogen | Fep. 22 | NP_001499.1. Peptide with sequence C-QSEAEPQSKSQ, from the internal region of the protein sequence according to

Peptide

C-QSEAEPQSKSQ Sequence

Shipping

Refrigerated Instructions

Predicted

Species

Human

Human

Gene ID

Product

2863

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

ELISA

Detection Limit

Antibody detection limit dilution 1:1000.

Western

Preliminary experiments in Human Brain (cerebellum) lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any Blot

results been reported with other antibodies/lysates?

Application Type

Pep-ELISA