

GOAT ANTI-GPR40 ANTIBODY

SKU: EB06761

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	GPCR40 FFA1R free fatty acid receptor 1 FFAR1 HGNC:4498 G protein-coupled receptor 40 GPR40
Accession ID	NP_005294.1
Blocking Peptide	EBP06761
Immunogen	Peptide with sequence CAARTQGGKSQK, from the C Terminus of the protein sequence according to NP_005294.1.
Peptide Sequence	CAARTQGGKSQK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	2864
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Approx 28kDa band observed in Human Breast lysates (calculated MW of 31.5kDa according to NP_005294.1). Recommended concentration: 0.01-0.03µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

