

GOAT ANTI-GPRC5B ANTIBODY

SKU: EB05406



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	G protein-coupled receptor, family C, group 1, member B retinoic acid responsive gene protein RAIG-2 RAIG2 G protein-coupled receptor, family C, group 5, member B GPRC5B
Blocking Peptide	EBP05406
Immunogen	Peptide with sequence C-TAPPSHTGRHLW, from the C Terminus of the protein sequence according to .
Peptide Sequence	C-TAPPSHTGRHLW
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Reactive Species	Human
Human Gene ID	51704
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 35-40kDa band seen in Human Pancreas lysate. Recommended for use at 2-4µg/ml.