



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-GPRIN3 ANTIBODY

SKU: EB10210

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias GPRIN3|GRIN3|GPRIN family member 3|G protein-regulated inducer of neurite outgrowth 3|FLJ42625

Names

Accession NP 938022.2 ID

Blocking

EBP10210 **Peptide**

Immunogen Peptide with sequence SNAQHTCKED, from the internal region of the protein sequence according to NP 938022.2.

Peptide

SNAQHTCKED Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted Human

Species

Human **Gene ID**

285513

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:1000.

Limit

Application Pep-ELISA **Type**