


GOAT ANTI-GRASP1 / GRIPAP1 ANTIBODY

SKU: EB06369



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	GRIPAP1 GRASP1 GRASP-1 KIAA1167 DKFZp434P0630 MPMGp800B12492Q3 GRIP1 associated protein 1 GRIP-associated protein 1 MGC126593 MGC126595
Accession ID	NP_064522.3; NP_997555.1
Blocking Peptide	EBP06369
Immunogen	Peptide with sequence LREQTTGLAAELQ-C, from the N Terminus of the protein sequence according to NP_064522.3; NP_997555.1.
Product Comments	This antibody is expected to recognize both reported isoforms (as represented by NP_064522.3; NP_997555.1)
Peptide Sequence	LREQTTGLAAELQ-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	56850
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Tonsil and Brain (Cortex). Recommended concentration: 3.75µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Preliminary experiments gave no signal but low background in Human Brain lysates at up to 1µg/ml.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

