

GOAT ANTI-GRB10 ANTIBODY

SKU: EB05151



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	KIAA0207 GRB-IR MEG1 IRBP RSS growth factor receptor-bound protein 10 GRB10
Accession ID	NP_005302
Blocking Peptide	EBP05151
Immunogen	Peptide with sequence KLKHH CIRVAL, from the C Terminus of the protein sequence according to NP_005302.
Peptide Sequence	KLKHH CIRVAL
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Reactive Species	Human
Human Gene ID	2887
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 65kDa band seen in HeLa cell lysates. Recommended for use at 4-6µg/ml.