

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-GRB10 ANTIBODY**

**SKU:** EB05151

## **SPECIFICATIONS**

**Formulation** Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

**Alias** KIAA0207|GRB-IR|MEG1|IRBP|RSS|growth factor receptor-bound protein 10|GRB10

**Names** 

**Accession** NP\_005302 ID

**Blocking** 

EBP05151 **Peptide** 

Immunogen Peptide with sequence KLKHHCIRVAL, from the C Terminus of the protein sequence according to NP\_005302.

**Peptide** 

KLKHHCIRVAL Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

**Predicted** 

Human, Mouse **Species** 

Reactive

Human

**Species** 

Human 2887 Gene ID

**Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

Western

Western blot: Approx 65kDa band seen in HeLa cell lysates. Recommended for use at 4-6µg/ml. **Blot**