

GOAT ANTI-GRB14 ANTIBODY

SKU: EB05080



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias	growth factor receptor-bound protein 14 GRB14
Names	
Usage Summary	has to fail as unexplainable 40 and 30kDa bands
Accession ID	NP_004481
Immunogen	Peptide with sequence KLKHYCARIAL, from the C Terminus of the protein sequence according to NP_004481.
Peptide Sequence	KLKHYCARIAL
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	2888
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Western blot: Approx 50-55kDa, 40kDa and 30kDa bands seen in HeLa cell lysates. Recommended for use at 0.1-0.5µg/ml.