

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-GRB14 ANTIBODY

SKU: EB05080

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

Alias growth factor receptor-bound protein 14|GRB14

Names

Usage **Summary**

has to fail as unexplainable 40 and 30kDa bands

Accession

NP 004481

Immunogen Peptide with sequence KLKHYCARIAL, from the C Terminus of the protein sequence according to NP_004481.

Peptide

KLKHYCARIAL Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human

Species

Human

Reactive **Species**

Human

2888 Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade **ELISA**

Antibody detection limit dilution 1:32000. **Detection**

Limit

Western Western blot: Approx 50-55kDa, 40kDa and 30kDa bands seen in HeLa cell lysates. Recommended for use at

Blot $0.1-0.5\mu g/ml$.