

GOAT ANTI-GRB7 (C TERMINUS) ANTIBODY

SKU: EB05079



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	growth factor receptor-bound protein 7 GRB7
Names	
Accession ID	NP_005301.2; NP_001025173.1
Immunogen	Peptide with sequence LLRHCCTRVAL, from the C Terminus of the protein sequence according to NP_005301.2; NP_001025173.1.
Product Comments	Both variants represent identical protein
Peptide Sequence	LLRHCCTRVAL
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	2886
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Western blot: Approx 50 kDa band seen in Hela and Jurkat Lysate. Recommended for use at 1-3 µg/ml. Please note that although the predicted size of GRB7 is approximately 59kDa according to NP_005301, literature suggests (for e.g., Pero et al., [J. Biol. Chem, Mar 2002; 277: 11918 - 11926]) observe GRB7 migrating at approx 50-55kDa on SDS-PAGE.
Application Type	Pep-ELISA, WB