


GOAT ANTI-GRB7 (N TERMINUS) ANTIBODY

SKU: EB05417



250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	OTTHUMP00000164352 growth factor receptor-bound protein 7 GRB7
Accession ID	NP_005301.2
Blocking Peptide	EBP05417
Immunogen	Peptide with sequence ELDLSPPHLSSSPE-C, from the N Terminus of the protein sequence according to NP_005301.2.
Product Comments	This antibody recognises the longer isoform of this protein in A431 lysate.
Peptide Sequence	ELDLSPPHLSSSPE-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	2886
Mouse Gene ID	14786
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 65kDa band observed in lysates of cell line A431 (calculated MW of 59.7kDa according to NP_005301.2). Recommended concentration: 1-3µg/ml. Primary incubation was 1 hour.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

