

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-GRB7 (N TERMINUS) ANTIBODY**

**SKU:** EB05417



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa



Telephone: (650) 697-3600

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / **Alias** 

OTTHUMP00000164352|growth factor receptor-bound protein 7|GRB7

Names

**Accession** 

NP\_005301.2

**Blocking** 

EBP05417 **Peptide** 

Peptide with sequence ELDLSPPHLSSSPE-C, from the N Terminus of the protein sequence according to **Immunogen** 

NP 005301.2.

**Product Comments** 

This antibody recognises the longer isoform of this protein in A431 lysate.

**Peptide** 

Sequence

**ELDLSPPHLSSSPE-C** 

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

**Predicted** 

Human **Species** 

Reactive **Species** 

Human

Human

2886 Gene ID

Mouse

14786 **Gene ID** 

**Product** 

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

FLISA

**Detection** Antibody detection limit dilution 1:128000.

Limit

Western Approx 65kDa band observed in lysates of cell line A431 (calculated MW of 59.7kDa according to

Blot NP\_005301.2). Recommended concentration: 1-3µg/ml. Primary incubation was 1 hour.

**Application** 

Pep-ELISA, WB Type



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GALLERY IMAGES**

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa