

GOAT ANTI-GRIM19 ANTIBODY

SKU: EB07096



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	CGI-39 protein B16.6, CDA016, CGI-39, GRIM-19 GRIM19 cell death-regulatory protein GRIM19 NADH dehydrogenase (ubiquinone) 1 alpha subcomplex
Accession ID	NP_057049.3
Blocking Peptide	EBP07096
Immunogen	Peptide with sequence C-RTTEEALHASH, from the C Terminus of the protein sequence according to NP_057049.3.
Peptide Sequence	C-RTTEEALHASH
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	51079
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Preliminary experiments in Human Heart, Skeletal Muscle and Placenta lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA