

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-GROWTH ARREST-SPECIFIC 6 ANTIBODY

SKU: EB05016

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias AXL stimulatory factor|AXLLG|AXSF|growth arrest-specific 6|GAS6

Names

Accession NP 000811

Immunogen Peptide with sequence DITAHSCPPVEPAAA, from the C Terminus of the protein sequence according to NP_000811.

Peptide DITAHSCPPVEPAAA

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Instructions Refrigerated

Predicted

Human **Species**

Human

2621 **Gene ID**

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade

Western

Western blot: No signal obtained yet but low background observed in Jurkat lysate at upto 3µg/ml. Blot