

GOAT ANTI-GROWTH ARREST-SPECIFIC 6 ANTIBODY

SKU: EB05016



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	AXL stimulatory factor AXLLG AXSF growth arrest-specific 6 GAS6
Accession ID	NP_000811
Immunogen	Peptide with sequence DITAHSCPPVEPAAA, from the C Terminus of the protein sequence according to NP_000811.
Peptide Sequence	DITAHSCPPVEPAAA
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	2621
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
Western Blot	Western blot: No signal obtained yet but low background observed in Jurkat lysate at upto 3µg/ml.